

DISSECTING MICROSCOPE



Glass Agencies are manufacturer of **Dissecting Microscope** Under **EROSE** Brand.

SPECIFICATION :

- Stand : Stable & Robudt with round base.
- Body : Equipped with joint revolving arm for moving the magnifying lens over the stage.
- Focusing : Corse motion corse movement by Rack and Pinion.
- Eyepieces : Two eye pieces of 10x and 20x.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

